



### Nasopharyngeal Swab



Breakpoint

### Oropharyngeal Swab




Breakpoint

## Advantages

Medline nylon flocked swabs provide optimum specimen collection because of their ability to collect more specimen at the collection site and to more rapidly release the entire sample specimen.


- Swab head is flocked nylon, and swab shaft is made of ABS (Acrylonitrile butadiene styrene).
- Non-toxic
- Individually wrapped and sterile




Cat.No.	Description	/Pk	/Cs
MD300235	Oropharyngeal Specimen Collection Swab, Individually Wrapped, Sterile	500	5000
MD300245	Nasopharyngeal Specimen Collection Swab, Individually Wrapped, Sterile	500	5000

## ▲ Technical Drawing

Nasopharyngeal Specimen Collection Swab



Oropharyngeal Specimen Collection Swab



## ▲ Operating Instruction


### Nasopharyngeal Specimen Collection Swab:

- With the patient seated, if possible, tilt the head back 70 degrees.
- Gently insert swab into the nostril along the septum floor of the nose extending straight back until the posterior nasopharynx is reached (distance from nostrils to external opening of ear).
- Rotate the swab several times while the swab is in contact with the nasopharyngeal wall.



### Oropharyngeal Specimen Collection Swab:

- With the patient seated, if possible, tilt the head back 70 degrees.
- Gently insert swab into back of the throat and tonsillar area.
- Rub the swab over both tonsillar pillars and the posterior oropharynx, avoid touching the tongue, teeth and gums.



# ▲ Technical Data

<i>Physical and chemical properties</i>	<b>Information on basic physical and chemical properties</b>	
	<b>General information</b>	<i>Plastic Shaft</i>
	Appearance	<i>Solid</i>
	Colour	<i>Milk White</i>
	Odour	<i>Characteristic</i>
	PH	<i>N/A</i>
	<b>Change in condition</b>	
	Melting point/ Melting range	<i>130°C</i>
	Boiling point/ Boiling range	<i>N/A</i>
	Flash point	<i>Not applicable</i>
	Explosive properties	<i>Product does not present an explosion hazard</i>
	Density	<i>1.03 ~ 1.07 / 25°C</i>
	Solubility in / Miscibility with water	<i>Insoluble</i>
<b>Toxicity</b>		
<i>No further relevant information available</i>		

<i>Physical and chemical properties</i>	<b>Information on basic physical and chemical properties</b>	
	<b>General information</b>	<i>Nylon Tip (Electrostatic Flocking)</i>
	Appearance	<i>Bulk Fibre</i>
	Colour	<i>White</i>
	Odour	<i>Characteristic</i>
	PH	<i>7</i>
	<b>Change in condition</b>	
	Melting point/ Melting range	<i>250°C</i>
	Boiling point/ Boiling range	<i>300°C</i>
	Flash point	<i>Not applicable</i>
	Explosive properties	<i>Product does not present an explosion hazard</i>
	Density	<i>1.03 ~ 1.07 / 25°C</i>
	Solubility in / Miscibility with water	<i>Insoluble</i>
<b>Toxicity</b>		
<i>No further relevant information available</i>		

<i>Physical and chemical properties</i>	<b>Information on basic physical and chemical properties</b>	
	<b>General information</b>	<i>PVC Flocking Glue</i>
	Appearance	<i>Liquid</i>
	Colour	<i>Light yellow</i>
	Odour	<i>Characteristic</i>
	PH	<i>7</i>
	<b>Change in condition</b>	
	Melting point/ Melting range	<i>5°C</i>
	Boiling point/ Boiling range	<i>121.2°C</i>
	Flash point	<i>Not applicable</i>
	Explosive properties	<i>Product does not present an explosion hazard</i>
	Density	<i>1.2</i>
	Solubility in / Miscibility with water	<i>Insoluble</i>
<b>Toxicity</b>		
LD50: 1535mg/kg :		
LC50: 50424 mg/m3		
<b>Composition</b>		
Cyclohexanone	2%	108-94-1
Xylene	5%	127-18-4



**medline**  
— s c i e n t i f i c —

**Head Office**

Medline Scientific Limited  
Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ  
United Kingdom

Tel: 00 44 1865 400321

Fax: 00 44 1865 400736

Email: [sales@medlinescientific.com](mailto:sales@medlinescientific.com)

[www.medlinescientific.com](http://www.medlinescientific.com)