

TECHNICAL DATA SHEET

Tube Deltaswab with Cary Blair media and standard flocked swab

code 304260KF

code	description	case quantity	case weight (kg)	case volume (m ³)
304260KF	2ml Cary Blair + standard flocked swab	6X100	6.30	0.056

Tube with media: Directive 98/79/CE. "In vitro" diagnostic medical devices.

Swab: Directive 93/42/CEE. Medical Devices.



INTENDED USE - DESCRIPTION

The kit includes:

- Tube with 2 ml Cary Blair media
- Standard flocked swab
- Flow-pack package

Sterile by radiation (R).

Especially suitable for the collection and transport of fecal samples.

Maintains viability of faecal pathogens for a minimum of 48 hours without overgrowth.

Designed to collect the sample directly from the rectum or of faeces. Designed for traditional and automated sample streaking.

Compatible with molecular diagnostic techniques by PCR and for direct extension on slides.

Tube dimensions

length	tube diameter	cap diameter
82.6 mm	13 mm	16.6 mm

Swab dimensions

Length: 152 mm

WRITTEN	VERIFIED	AUTHORISED
Joan Ibáñez Marketing Department	Anna Mir Technical, Development, Quality and Environment Director	Mónica Torras Sales & Marketing Director

TECHNICAL DATA SHEET

Tube Deltaswab with Cary Blair media and standard flocked swab

code 304260KF


FEATURES AND ADVANTAGES OF DELTASWAB SYSTEM

- Compatible with the new inoculation automation and sample streaking equipment.
- Compatible with molecular diagnostic techniques.
- Facilitates the collection, transport and subsequent treatment of microbiological samples.
- Greater reliability for sample recovery: high capacity for the absorption and elution of the sample from the flocked swab deltaswab.
- Greater feasibility of the sample as it is totally suspended in the medium.
- Allows different homogeneous inocula to be obtained from the same sample.
- It adapts to any working protocol.
- Minimises any possible cross-contamination: less handling and maximum water tightness.
- More comfort for the user: conical base which facilitates stirring and greater stability of the tube with skirt base.
- More comfort for the patient, due to the softer swab covering.
- Standardises the receipt of samples at the microbiology labs.
- Facilitates manual streaking.
- Allows direct extension on the glass slide as the medium does not contain agar.
- Ensures proper sample transport and storage both at room temperatures (20 °C – 25 °C) and at refrigeration temperatures (4 °C – 8 °C).
- Maximum reliability and convenience of the breakpoint.

FEATURES AND ADVANTAGES OF FLOCKED SWABS

- Support made of polystyrene --> Optimum safety.
- Head swab made of Flocked fibers --> increased absorption and elution of the sample on the media.
- Maximum adaptability: wide range
- Maximum comfort: anatomical and ergonomic design.

INSTRUCTIONS



1. Open flow-pack using both hands pulling opposite sides.
2. Collect the sample with the swab.
3. Open the tube with the other hand and place the swab inside so that it is covered with the media. Align the breakpoint of the swab with the top of the tube lightly pressing the swab downwards.
4. Break the swab by the breakpoint by supporting it on the inner edge of the tube.
5. Discard leftover piece of stick, screw cap tightly and shake the sample to elute it into the media.

Send the tube with the sample to the laboratory for analysis

WRITTEN	VERIFIED	AUTHORISED
Joan Ibáñez Marketing Department	Anna Mir Technical, Development, Quality and Environment Director	Mónica Torras Sales & Marketing Director

TECHNICAL DATA SHEET


Tube Deltaswab with Cary Blair media and standard flocked swab

code 304260KF

TUBE AND SWAB KIT

code	description
304260KF	2ml Cary Blair + standard flocked swab

LABEL DESIGN




WRITTEN	VERIFIED	AUTHORISED
Joan Ibáñez Marketing Department	Anna Mir Technical, Development, Quality and Environment Director	Mónica Torras Sales & Marketing Director

TECHNICAL DATA SHEET

Tube Deltaswab with Cary Blair media and standard flocked swab

code 304260KF

FLOW-PACK DESIGN



WRITTEN	VERIFIED	AUTHORISED
Joan Ibáñez Marketing Department	Anna Mir Technical, Development, Quality and Environment Director	Mónica Torras Sales & Marketing Director