

## White 6.7mm Phthalate Free PU ESD Medical Gas Hose (Oxygen)

Product Code: I08601PF

The White 6.7mm Phthalate Free PU ESD Medical Gas Hose (Oxygen) is designed to reduce the risk of explosions from modern anaesthetic gases.

All the materials used are USP VI compliant.



### Specifications

<b>Colour</b>	White
<b>Coil Length (m)</b>	30
<b>Internal Diameter</b>	6.7
<b>Gas Type</b>	Oxygen

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

Tel: (+44) 01865 400321

Email: [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

Website: [www.medlinescientific.com](http://www.medlinescientific.com)