

## Yellow 6.7mm Phthalate Free PU ESD Gas Medical Hose (Vacuum)

Product Code: I08603PF

The Yellow 6.7mm Phthalate Free PU ESD Gas Medical Hose (Vacuum) is designed to reduce the risk of explosions from modern anaesthetic gases.

All the materials used are USP VI compliant.



### Specifications

<b>Colour</b>	Yellow
<b>Coil Length (m)</b>	30
<b>Internal Diameter (mm)</b>	6.7
<b>Gas Type</b>	Vacuum

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)