

Black/White 6.7mm Striped Phthalate Free PU ESD Gas Medical Hose (Air)

Product Code: I086I2PF

The Black/White 6.7mm Striped Phthalate Free PU ESD Medical Hose (Air) is designed to reduce the risk of explosions from modern anaesthetic gases.

All the materials used are USP VI compliant.



Specifications

Colour	Black/White (striped)
Coil Length (m)	30
Internal Diameter (mm)	6.7
Gas Type	Air

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com