

Yellow 12.7mm Phthalate Free PU ESD Gas Medical Hose (Vacuum)

Product Code: I08643PF

The Yellow 12.7mm Phthalate Free PU ESD Gas Medical Hose (Vacuum) is designed to reduce the risk of explosions from modern anaesthetic gases.

All the materials used are USP VI compliant.



Specifications

Colour	Yellow
Coil Length (m)	30
Internal Diameter (mm)	12.7
Gas Type	Vacuum

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com