

## **CETI Colt Trinocular Compound Microscope (4x, 10x, 40x, 100x Achromatic)**

**Product Code: 1100.5103**

**Brand: CETI**

Designed to provide students with a high performance, CETI's Colt Trinocular Compound Microscope (4x, 10x, 40x, 100x Achromatic) is an excellent choice for colleges/universities.

Its 10x/16mm wide-field eyepiece enables distinctive points on slides to be seen quickly and efficiently. With its achromatic objectives, this educational microscope is guaranteed to produce sharp, colour-corrected images.



## Specifications:

<b>Optical Head</b>	Trinocular, inclined at 30° and 360° (rotatable)
<b>Eye Piece</b>	10x/16mm wide-field
<b>Nosepiece</b>	Quadruple, revolving with click stop
<b>Objectives</b>	4x, 10x, 40x, 100x Achromatic
<b>Condenser</b>	Abbe condenser with iris diaphragm
<b>Stage</b>	Size 120 x 120mm
	<b>Moving Range:</b> 60 x 30mm
<b>Focusing Knobs</b>	Coaxial coarse and fine
<b>Mains Power</b>	220 - 240V/50 - 60Hz (Euro connector)
<b>Converter</b>	Built-in low voltage
<b>Bulb</b>	Halogen lamp 6V/15W/G4
<b>Supplied With</b>	One spare bulb, two spare fuses, dust cover, and instruction manual
<b>Packing</b>	Two Styrofoam shelves in cardboard box

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)  
**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)