

CETI Colt Monocular Compound Microscope (4x, 10x, 40x, 60x Achromatic)

Product Code: I100.5160M

Brand: CETI

The affordable CETI Colt Monocular Compound Microscope (4x, 10x, 40x, 60x Achromatic) offers undergraduate and postgraduate students flawless visibility of samples and brilliant optical quality for various laboratory applications.

This educational microscope comes with an Abbe condenser with iris diaphragm, along with coaxial coarse and fine focusing knobs. The magnification of this compound microscopes ranges from 4x to 60x.



Specifications:

Optical Head	Monocular, inclined at 30° and 360° (rotatable)
Eye Piece	10x/16mm wide-field
Nosepiece	Quadruple, revolving with click stop
Objectives	4x, 10x, 40x, 60x Achromatic
Condenser	Abbe condenser with iris diaphragm
Stage	Size 120 x 120mm
	Moving Range: 60 x 30mm
Focusing Knobs	Coaxial coarse and fine
Mains Power	220 - 240V/50 - 60Hz (Euro connector)
Converter	Built-in low voltage
Bulb	Halogen lamp 6V/15W/G4
Supplied With	One spare bulb, two spare fuses, dust cover, and instruction manual
Packing	Two Styrofoam shelves in cardboard box

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com
Website: www.medlinescientific.com