

## **Ceti Akropol Monocular Polarising Compound Microscope (4x, 10x, 40x, 60x Achromatic)**

**Product Code: 1111.2441**

**Brand: Ceti**

The Ceti Akropol Monocular Polarising Compound Microscope (4x, 10x, 40x, 60x Achromatic) is designed to introduce polarisation techniques, whilst allowing for the measurement of the degree of rotation that samples induce.

It is supplied with one spare bulb, spare fuse, immersion oil (10ml), dust cover, and instruction manual.



## Specifications

<b>Optical Head</b>	Monocular
<b>Eyepiece</b>	High eye point 10x/16mm wide-field
<b>Objectives</b>	4x, 10x, 40x, 60x Achromatic
<b>Condenser</b>	Brightfield Abbe - NA 1.25
<b>Stage</b>	Round centrable and rotating stage 360°
<b>Focusing</b>	Coaxial coarse and fine
<b>Illumination (variable)</b>	Halogen 20W

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)