

Ceti Max III Trinocular Compound Microscope (4x, 10x, 40x, 100x Infinity Corrected Plan Achromatic)

Product Code: I203.3006

Brand: Ceti

Ceti's Max III Trinocular Compound Microscope (4x, 10x, 40x, 100x Infinity Corrected Plan Achromatic) is professionally designed for laboratory applications in teaching laboratories, offering precise high-quality images with its fully-coated optical system.

This microscope comes with an inclined 30° and 360° rotatable trinocular head and a built-in mechanical stage. Four high clarity Infinity Corrected Plan Achromatic objectives are supplied as standard.



INTEGRAL CARRYING HANDLE

Allows for one handed transportation of the microscope.



OBJECTIVES

A wide variety of high clarity objectives are available with Infinity Corrected Plan Achromatic optics.

All objectives are colour coded to ease identification when in use.



PRECONDENSER

Fixed lens, which is removable to accept a mirror (for locations without power) and to allow for lamp replacement.



REVERSED NOSEPIECE & FULL MECHANICAL STAGE

Quadruple reversed, so providing easy uninterrupted access to the stage when changing over slides. Has full X/Y movement and spring loaded slide holder.

Specifications

Optical Head	Trinocular with photo port, inclined at 30°
	Interpupillary distance adjustment (Siedentopf type): 47 to 78mm
	Dioptric adjustment on one eyepiece
Eyepieces	10x/18mm wide-field (pointer in one eyepiece)
Nosepiece	Quadruple reversed, revolving with click stop
Objectives	4x, 10x, 40x, 100x Infinity Corrected Plan Achromatic
Condenser	Brightfield Abbe NA 1.25 - precentred, focusable
	Iris diaphragm and filter holder
Stage	160 x 140mm with built-in mechanical stage and slide holder
	75 x 50mm movement, Vernier scale
Focusing Knobs	Coaxial coarse and fine with stop
Mains Power	220 - 240V / 50 - 60Hz
Converter	Built-in low voltage (3.2 - 3.5V electric)
Illumination	LED 3W with intensity control (Halogen also available)
Supplied With	Immersion oil (5ml), blue filter, dust cover, power cord with unique 'self-store' system, and instruction manual

Accessories

CODE	DESCRIPTION	CODE	DESCRIPTION
	Accessories		
1202.4016	Eye piece single, 16x / 18mm wide-field	1202.1401	Aluminium case
1202.4017	Eye piece pair, 16x / 18mm wide-field	1203.2202	Mirror, plan concave, with holder
1202.1005	Eye piece pair 10x/18mm wide-field	1202.4888	C mount adapter
2249.1001	Object micrometer with scale 1mm 0.01mm used with micrometer eyepiece	2280.0101	Immersion oil 5 ml
1202.4007	Eye piece micrometer 10x/18mm with cross line scale 0-18mm 0.1mm (single)	2480.2075	Cleaning set lens cleaning tissues (20 pcs) + solution (75ml) + brush
1202.4020	Objective achromatic 20x / 0.40 / 160 / 0.17		Spare Parts
1202.4066	Objective achromatic 60x / 0.85 / 160 / 0.17 (dry) spring loaded	2260.4000	Bulb halogen 6V/20W/G4 Each
1202.5507	Polarising set, analyser + polariser	2260.5000	LED 5W replacement.
1202.4603	Darkfield condenser NA 0.9	1202.4004	Objective achromatic 4x / 0.10 / 160 / 0.17
1202.4703	Darkfield condenser oil NA1.25-1.36	1202.4010	Objective achromatic 10x / 0.25 / 160 / 0.17
1202.4999	Turret phase contrast kit with Ach 10x, 40x 100x phase objectives	1202.4040	Objective achromatic 40x / 0.65 / 160 / 0.17 spring loaded
1202.4777	Simple phase contrast kit with Ach 10x 40x phase objectives and separate condenser	1202.4100	Objective achromatic 100x / 1.25 / 160 / 0.17 (oil) spring loaded
2460.2002	Heating stage 150x130cm Amb to 50° C including controller	1202.4601	Condenser Abbe (NA 1.25) with iris diaphragm + filter holder
		1202.9212	Filter blue Ø 32mm
		2299.1524	Fuse 5x 20/500 MA (10 pcs)

Unit 3, Tower Business Park
 Warpsgrove Lane
 Chalgrove
 Oxfordshire
 OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com
 Website: www.medlinescientific.com