

## Deltalab Sterilisation Indicator Tape (Ethylene Oxide)

**Product Code:** I9I050

**Brand:** Deltalab

Deltalab's Sterilisation Indicator Tape (Ethylene Oxide) is self-adhesive and made of semi-creped paper. It also sticks on different surfaces such as plastic, metal, fabrics, paper, cardboard, etc. Will remain affixed for 18 months.

### Features:

- Thickness range from 0.15 to 0.17mm.
- Indicator for sterilisation by ethylene oxide gas. The word GAS printed in brown changes to green once it has been sterilised as per the following specifications: exposition to 500mg/l EO and 60%RH for minimum 30 minutes at 30°C, and exposition to 12% EO and 60%RH for minimum 48 minutes at 54°C.

### Specifications:

|                       |                    |
|-----------------------|--------------------|
| Type of Sterilisation | Ethylene Oxide Gas |
| Length (mm)           | 50                 |
| Width (mm)            | 19                 |
| Sterility             | Non-Sterile        |
| Quantity              | 1                  |



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)