

## Deltalab Sterilisation Indicator Tape (Dry Heat - Poupinel)

**Product Code:** I91051

**Brand:** Deltalab

Deltalab's Sterilisation Indicator Tape (Dry Heat - Poupinel) is self-adhesive and made of semi-creped paper. It also sticks on different surfaces such as plastic, metal, fabrics, paper, cardboard, etc. Will remain affixed for 18 months.

### Features:

- Thickness range from 0.15 to 0.17mm.
- Indicator for sterilisation with poupinel method (dry heat). The green strips turn to brown by the effect of time and temperature as per the following specification: exposition a minimum 20 minutes at 160°C or more in dry air.

### Specifications:

<b>Type of Sterilisation</b>	Dry Heat (Poupinel)
<b>Length (mm)</b>	50
<b>Width (mm)</b>	19
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	1

