

Deltalab Sterilisation Indicator Tape (Autoclave)

Product Code: I9I052

Brand: Deltalab

Deltalab's Sterilisation Indicator Tape (Autoclave) is self-adhesive and made of semi-creped paper. It also sticks on different surfaces such as plastic, metal, fabrics, paper, cardboard, etc. Will remain affixed for 18 months.

Features:

- Thickness range from 0.15 to 0.17mm.
- Indicator for sterilisation with steam. The green strips turn to dark brown by the effect of time, temperature, and saturated steam as per the following specifications: exposition at 121°C for 10 minutes and at 134°C for 2 minutes.

Specifications:

Type of Sterilisation	Autoclave
Length (mm)	50
Width (mm)	19
Sterility	Non-Sterile
Quantity	1



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com