

Deltalab Measuring Scoops (10ml)

Product Code: I91057

Brand: Deltalab

Made with materials suitable for alimentary use, Deltalab's Measuring Scoops (10ml) are autoclavable and have a flat base for stability, along with a flat top for levelling powders.

Specifications:

Alimentary Use?	Yes
Capacity (ml)	10
Length (mm)	100
Material	Polypropylene
Sterility	Non-Sterile
Quantity	12



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com