

## Dentalab Measuring Scoops (25ml)

**Product Code:** I91058

**Brand:** Dentalab

Made with materials suitable for alimentary use, Dentalab's 25ml Measuring Scoops are autoclavable and have a flat base for stability, along with a flat top for levelling powders.

### Specifications:

<b>Alimentary Use?</b>	Yes
<b>Capacity (ml)</b>	25
<b>Length (mm)</b>	135
<b>Material</b>	Polypropylene
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	12



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)