

Deltalab Coplin Type Staining Jar

Product Code: I91087

Brand: Deltalab

Deltalab's Coplin Type Staining Jar is made of opaque white, unbreakable and autoclavable polypropylene. It is designed to stain 5 standard microscope slides.

It features a rectangular container with a rounded base and screw ribbed flat cap with an inner ring for preventing leakage.

Used as a storage system, it can hold up to 10 slides (2 per groove).

Specifications:

Base Diameter (mm)	60
Height of Jar (mm)	110
Height of Jar with Cap (mm)	114
Material	Polypropylene
Case Quantity	12

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com