

## Deltalab White Diamond Weighing Dishes (41 x 41 x 8mm)

**Product Code:** 1911103


**Brand:** Deltalab

Deltalab's White Diamond Weighing Dishes (41 x 41 x 8mm) are made of translucent polystyrene and can be easily bent into pouring spouts and enable non-spill transfer. Biologically inert and contaminant-free, these weighing dishes are resistant to diluted acids, aqueous solutions, alcohols, and bases. They are also suitable for microwaves and non-antistatic.

They will withstand temperatures from 80°C to -10°C.

### Specifications:

<b>Dimensions (mm)</b>	41 x 41 x 8
<b>Weight (g)</b>	0.34
<b>Material</b>	Polystyrene
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	1000



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)