

Deltalab White Diamond Weighing Dishes (140 x 140 x 22mm)

Product Code: 1911105


Brand: Deltalab

Deltalab's White Square Weighing Dishes (140 x 140 x 22mm) are made of translucent polystyrene and can be easily bent into pouring spouts and enable non-spill transfer. Biologically inert and contaminant-free, these weighing dishes are resistant to diluted acids, aqueous solutions, alcohols, and bases. They are also suitable for microwaves and non-antistatic.

They will withstand temperatures from 80°C to -10°C.

Specifications:

Dimensions (mm)	140 x 140 x 22
Weight (g)	3.27
Material	Polystyrene
Sterility	Non-Sterile
Quantity	250



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com