

Deltalab Analytical Funnels (45mm)

Product Code: I9147

Brand: Deltalab

Made of autoclavable polypropylene, Deltalab's Analytical Funnels (45mm) feature a cone angle at 60°C.

Specifications:

Diameter (mm)	45
Stem Diameter (mm)	5
Stem Length (mm)	43
Total Length (mm)	79
Material	Polypropylene
Quantity	20

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com