

Deltalab Yellow Wide Mouth Wash Bottles (250ml)

Product Code: I91637.06

Brand: Deltalab

Deltalab's 250ml Red Capped Wide Mouth Wash Bottles are translucent low density polyethylene bottles with a screw cap and delivery tube. Suitable for a variety of liquid storage and sampling applications, the tops of these bottles provide easy identification. The wide mouth is used to provide fast and easy filling.

Features:

- Feature a one piece cap and delivery tube assembly, thus eliminating air or liquid absorption through the cap.
- Includes a wide neck for safe and easy filling.
- Ribbed caps for a better handling while opening and closing.
- Graduated in 25ml increments.

Specifications:

Capacity (ml)	250
Cap Colour	Red
Internal Neck (Ø mm)	38
Base (Ø mm)	60
Height (mm)	140
Material	Polyethylene
Sterility	Non-Sterile
Quantity	10



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com