

## Deltalab Wide Neck Bottles (10 Litres)

**Product Code:** 191642

**Brand:** Deltalab

Made of high density polyethylene, Deltalab's Wide Neck Bottles (10 Litres) feature an insert plug for a positive leakproof seal and one handle for easy carrying.

There are also tagging points on the neck and cap for a tamper proof closure or for labelling.

They are designed for both liquids and powders, and are easy to fill and pour.

### Specifications:

<b>Capacity (litres)</b>	10
<b>Neck Diameter (mm)</b>	120
<b>Body Diameter (mm)</b>	210
<b>Material</b>	Polyethylene
<b>Quantity</b>	1



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)