

Deltalab Drying Rack

Product Code: I92I3

Brand: Deltalab

Deltalab's Drying Rack is a wall mountable high impact white polystyrene drying rack for general labware. It consists of a single mould 4mm plate with 72 peg sockets.

The push-in pegs can be placed in any configuration on the rack and are easily removable to accommodate odd shaped items.

The peg sockets have closed ends to avoid dripping and eliminate the risk of biological contamination.

The rack is suitable for items with a mother diameter larger than 15mm.

A drip channel collects dripping that drains via a drain tube provided. It includes 11 additional peg sockets (Ø 6mm) for smaller parts.

The rack comes complete with a mounting kit and fixing template.

Specifications:

Dimensions (mm)	450 x 630 x 110
Case Quantity	1



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com