

500ml Graduated narrow neck bottles

Product Code: I9324

Brand: Deltalab

Made of translucent polyethylene.
Designed for the storage of acids and bases.
Include a ribbed screw cap for a better handling, and an insert plug.
Feature sloping shoulders for easy pouring.
Tagging points on neck and cap for a tamper proof closure and for labeling.
Bottles have the molded graduations.

Supplied capped.

Specifications:

Capacity (ml)	500
Neck Ø (mm)	23
Bottom Ø (mm)	74
Height (mm)	157
Sterility	Non-Sterile
Quantity	6 x 25



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com