

Deltalab Buchner Funnel (810ml)

Product Code: I9442

Brand: Deltalab

Made of autoclavable polypropylene, Deltalab's Buchner Funnel (810ml) is an excellent and reliable item for carrying out filtration applications in laboratories.

It is unbreakable and resistant to corrosion and heat. It can also be easily dismantled for cleaning.



Specifications:

Capacity (ml)	810
Filter Diameter (mm)	110
Upper Diameter (mm)	30.1
Bottom Diameter (mm)	22.6
Holes (mm)	2
Material	Polypropylene
Case Quantity	1

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com