

## Deltalab Buchner Funnel (810ml)

**Product Code:** I9442

**Brand:** Deltalab

Made of autoclavable polypropylene, Deltalab's Buchner Funnel (810ml) is an excellent and reliable item for carrying out filtration applications in laboratories.

It is unbreakable and resistant to corrosion and heat. It can also be easily dismantled for cleaning.

### Specifications:

<b>Capacity (ml)</b>	810
<b>Filter Diameter (mm)</b>	110
<b>Upper Diameter (mm)</b>	30.1
<b>Bottom Diameter (mm)</b>	22.6
<b>Holes (mm)</b>	2
<b>Material</b>	Polypropylene
<b>Case Quantity</b>	1

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)