


Deltalab Automatic Slide Dispenser

Product Code: I9540

Brand: Deltalab

Made of ABS and PS, Deltalab's Automatic Slide Dispenser can hold up to 50 standard slides. Slides can be dispensed individually by rotating the knurled knob on either side of the moulded base.



Specifications:

Description	For 50 slides (26 x 76mm)
Case Quantity	1

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com