

## Deltalab Automatic Slide Dispenser

**Product Code:** I9540

**Brand:** Deltalab

Made of ABS and PS, Deltalab's Automatic Slide Dispenser can hold up to 50 standard slides. Slides can be dispensed individually by rotating the knurled knob on either side of the moulded base.

### Specifications:

<b>Description</b>	For 50 slides (26 x 76mm)
<b>Case Quantity</b>	1

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)