

## Universal test tube racks


**Product Code:** I9564.06

**Brand:** Deltalab

Made of autoclavable polypropylene.  
Specially designed to ensure that all test tubes will be incubated at the same temperature.  
An alpha-numeric grid reference on top tier of the rack allows an easy sample location.  
Racks are easily assembled, and are stackable when unloaded.  
Each rack is supplied dismantled.  
Easy to assemble.

### Specifications:

<b>Places</b>	90 (6x15)
<b>Test tube Ø (mm)</b>	23
<b>Dimensions (mm)</b>	80 x 110 x 45
<b>Colour</b>	Yellow
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	5



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)