

Deltalab Pipette Rack

Product Code: I9973

Brand: Deltalab

Deltalab's Pipette Rack is designed for orderly, space-saving pipette storage. This laboratory rack holds 50 pipettes in a vertical 5 x 10 arrangement, ensuring that the pipettes will not get in contact with the table surface.

The end plates have handles for easy carrying.

Specifications:

Material	Autoclavable Polypropylene
Hole Diameter (mm)	16
Dimensions (mm)	213 x 114 x 222
Case Quantity	6

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com