

## Deltalab Ethylene Oxide Sterilised Graduated Pipettes (5ml) - Peel-Packed

**Product Code:** 200004

**Brand:** Deltalab

Deltalab's Ethylene Oxide Sterilised Graduated Pipettes (5ml) - Peel-Packed) are made of transparent, inert, and low density polyethylene. Designed for fast liquid transfer, these pasteur pipettes each have a built-in bulb and metric graduations that allow liquids and suspensions to be handled easily and efficiently.



### Specifications:

<b>Volumetrical Tolerance</b>	±10%
<b>Quantity</b>	6 x 500 (individually wrapped, peel-packs)
<b>Length (Ø mm)</b>	154
<b>Stem (Ø mm)</b>	5
<b>Graduated Until (ml)</b>	1.5
<b>Total Capacity (ml)</b>	5
<b>Bulb Draw (ml)</b>	3.1
<b>Material</b>	Polyethylene
<b>Graduations (ml)</b>	0.25
<b>Sterility</b>	Ethylene Oxide

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)