


Deltalab Ethylene Oxide Sterilised Graduated Pipettes (5ml) - Peel-Packed

Product Code: 200004

Brand: Deltalab

Deltalab's Ethylene Oxide Sterilised Graduated Pipettes (5ml) - Peel-Packed) are made of transparent, inert, and low density polyethylene. Designed for fast liquid transfer, these pasteur pipettes each have a built-in bulb and metric graduations that allow liquids and suspensions to be handled easily and efficiently.



Specifications:

Volumetrical Tolerance	±10%
Quantity	6 x 500 (individually wrapped, peel-packs)
Length (Ø mm)	154
Stem (Ø mm)	5
Graduated Until (ml)	1
Total Capacity (ml)	5
Bulb Draw (ml)	3.1
Material	Polyethylene
Graduations (ml)	0.25
Sterility	Ethylene Oxide

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com