

Deltalab Yellow Eppendorf Pipette Tip Rack (5 - 200µl)

Product Code: 200009R

Brand: Deltalab

Deltalab's Yellow Eppendorf Pipette Tip Rack (5 - 200µl) features a transparent lid and blue base, which can hold up to 96 yellow pipette tips.

It is made of autoclavable polypropylene up to 121°C.

Specifications:

Description	96 tips rack
Sterility	Sterile
Quantity	10



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com