

## Deltalab Natural Colour Eppendorf Pipette Tip Rack (2 - 20µl)

**Product Code:** 200010FR

**Brand:** Deltalab

Deltalab's Natural Colour Eppendorf Pipette Tip Rack (2 - 20µl) has a transparent hinged lid and base. It can hold up to 96 pipette tips.

It is made of autoclavable polypropylene up to 121°C.



### Specifications:

<b>Description</b>	96 tips rack filter
<b>Material</b>	Polypropylene
<b>Sterility</b>	Sterile
<b>Quantity</b>	10

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)