

## Deltalab Yellow Gilson Pipette Tip Rack (2 - 200µl)

**Product Code:** 200016R

**Brand:** Deltalab

Deltalab's Yellow Gilson Pipette Tip Rack (2 - 200µl) features a transparent lid and blue base and it can hold up to 96 yellow pipette tips.

The rack is made of autoclavable (up to 121°C) polypropylene.

### Specifications:

<b>Description</b>	96 tips rack
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	10



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)