

## Deltalab Natural Colour Universal Tips (1 - 200µl)


**Product Code:** 200026

**Brand:** Deltalab

Deltalab's Natural Colour Universal Tip (1 - 200µl) are moulded accurately to provide a fine fit and performance for a range of different pipettors. Each tip is manufactured in high-quality polypropylene.

### Specifications:

<b>Description</b>	Tips in bulk
<b>Length (mm)</b>	91
<b>Material</b>	Polypropylene
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	5 x 1000



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)