


## Deltalab Ethylene Oxide Sterilised Graduated Pipettes (7ml) - Packs of 20

**Product Code:** 200062R

**Brand:** Deltalab

Deltalab's Ethylene Oxide Sterilised Graduated Pipettes (7ml) Packs of 20 are made of transparent, inert, and low density polyethylene. Designed for fast liquid transfer, these pasteur pipettes each have a built-in bulb and metric graduations that allow liquids and suspensions to be handled easily and efficiently.



### Specifications:

<b>Volumetrical Tolerance</b>	±10%
<b>Quantity</b>	10 x 500 (packs of 20)
<b>Length (Ø mm)</b>	154
<b>Stem (Ø mm)</b>	7.8
<b>Graduated Until (ml)</b>	3
<b>Total Capacity (ml)</b>	7
<b>Bulb Draw (ml)</b>	3.2
<b>Material</b>	Polyethylene
<b>Graduations (ml)</b>	0.50
<b>Sterility</b>	Ethylene Oxide

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)