

Deltalab Natural Colour Pipette Tip Rack (100 - 1000µl)

Product Code: 200070NR

Brand: Deltalab

Deltalab's Natural Colour Pipette Tip Rack (100 - 1000µl) has a transparent lid and blue base. It can hold up to 96 Eppendorf®, Gilson®, Biohit®, DSG®, Elkay®, Finnipipette®, Genex®, Jencons®, Nichiryo®, and Menarini® tips.

Made of polypropylene, this rack is autoclavable up to 121°C.

Specifications:

Description	96 tips rack
Material	Polypropylene
Sterility	Sterile
Quantity	10



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com