

Deltalab Natural Colour Eppendorf Pipette Tip Rack (10 - 300µl)

Product Code: 200077FR

Brand: Deltalab

Deltalab's Natural Colour Eppendorf Pipette Tip Rack (10 - 300µl) has a hinged transparent lid and base. It can hold up to 96 tips.

This rack is made of autoclavable (up to 121°C) polypropylene.

Specifications:

Description	96 tips rack filter
Material	Polypropylene
Sterility	Sterile
Quantity	10



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com