

Deltalab Coulter Counter Cups Without Lids (20ml)

Product Code: 200102

Brand: Deltalab

Used for cell counting, Deltalab's Coulter Counter Cups Without Lids (20ml) are single use polystyrene cups suitable for any Coulter.

Specifications:

Volume (ml)	20
Dimensions (Ø x h)	30 x 56mm
Description	Coulter counter cup without lid
Material	Polystyrene
Case Quantity	1400



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com