

## 90mm Triple Vent Petri Dish

**Product Code:** 200200.5

**Brand:** Medline

Made of Virgin polystyrene, Medline Petri Dishes Internal Crosses (90 x 15.8mm) come with three vents and are aseptically sterile. Suitable for automatic filling machines, these stackable petri dishes can be used in microbiological culturing on solid media



### Specifications:

<b>Dimensions (mm)</b>	90 x 15
<b>Material</b>	Polystyrene
<b>Sterility</b>	Sterile
<b>Quantity</b>	500

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)