

Deltalab Petri Dishes, Three Compartments (90 x 14mm)

Product Code: 200203

Brand: Deltalab

Made of polystyrene and aseptic, Deltalab's Petri Dishes, Three Compartments (90 x 14mm) are supplied in groups of 15 units, packaged in heat sealed bags. They are suitable for automatic filling machines.

Specifications:

Dimensions	Ø 90 x 14mm
Material	Polystyrene
Case Quantity	25 x 20



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com