

## Deltalab Sterile Contact Plates

**Product Code:** 200208

**Brand:** Deltalab

Deltalab's Sterile Contact Plates can be used in the pharmaceutical industry and hospitals for determining bacterial contamination of surfaces such as the skin, operating tables, refrigerated gondolas, and work surfaces.

Manufactured in transparent polystyrene, their moulded grid at the bottom makes it possible to identify the detected contamination per cm<sup>2</sup> and it facilitates the counting of colonies.

The shape of the lid makes stacking perfectly stable during transport and incubation, and saves space on work surfaces.

### Specifications:

<b>Dimensions</b>	65.7 x 14.7mm
<b>Description</b>	Roadac contact plate
<b>Sterility</b>	Non-Sterile
<b>Material</b>	Polystyrene
<b>Case Quantity</b>	25 x 20

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)