

Deltalab Sterile Contact Plates

Product Code: 200208

Brand: Deltalab

Deltalab's Sterile Contact Plates can be used in the pharmaceutical industry and hospitals for determining bacterial contamination of surfaces such as the skin, operating tables, refrigerated gondolas, and work surfaces.

Manufactured in transparent polystyrene, their moulded grid at the bottom makes it possible to identify the detected contamination per cm² and it facilitates the counting of colonies.

The shape of the lid makes stacking perfectly stable during transport and incubation, and saves space on work surfaces.



Specifications:

Dimensions	65.7 x 14.7mm
Description	Roadac contact plate
Sterility	Non-Sterile
Material	Polystyrene
Case Quantity	25 x 20

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com