

## Deltalab Autoclave Bags (30 x 60cm)

**Product Code:** 200310

**Brand:** Deltalab

Deltalab's Autoclave Bags (30 x 60cm) are printed with a biohazard graphic symbol and instructions of use in five languages. Do not close the bags hermetically when placing them in the autoclave at 121°C.

### Specifications:

<b>Dimensions (cm)</b>	30 x 60
<b>Material</b>	Polypropylene
<b>Approx. Capacity (L)</b>	12
<b>Thickness (microns)</b>	37.5
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	500



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)