

Deltalab Autoclave Bags (30 x 60cm)

Product Code: 200310

Brand: Deltalab

Deltalab's Autoclave Bags (30 x 60cm) are printed with a biohazard graphic symbol and instructions of use in five languages. Do not close the bags hermetically when placing them in the autoclave at 121°C.

Specifications:

Dimensions (cm)	30 x 60
Material	Polypropylene
Approx. Capacity (L)	12
Thickness (microns)	37.5
Sterility	Non-Sterile
Quantity	500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com