

## Deltalab Autoclave Bags (40 x 75cm)

**Product Code:** 200311

**Brand:** Deltalab

Deltalab's Autoclave Bags (40 x 75cm) are printed with a biohazard graphic symbol and instructions of use in five languages. Do not close the bags hermetically when placing them in the autoclave at 121°C.

### Specifications:

|                             |               |
|-----------------------------|---------------|
| <b>Dimensions (cm)</b>      | 40 x 75       |
| <b>Material</b>             | Polypropylene |
| <b>Approx. Capacity (L)</b> | 24            |
| <b>Thickness (microns)</b>  | 37            |
| <b>Sterility</b>            | Non-Sterile   |
| <b>Quantity</b>             | 300           |



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)