

Deltalab Autoclave Bags (60 x 75cm)

Product Code: 200312

Brand: Deltalab

Deltalab's Autoclave Bags (60 x 75cm) are printed with a biohazard graphic symbol and instructions of use in five languages. Do not close the bags hermetically when placing them in the autoclave at 121°C.

Specifications:

Dimensions (cm)	60 x 75
Material	Polypropylene
Approx. Capacity (L)	49
Thickness (microns)	37
Sterility	Non-Sterile
Quantity	200

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com