

Deltalab Autoclave Resistant Red Bags (48 x 58cm)

Product Code: 200320

Brand: Deltalab

Deltalab's Autoclave Resistant Red Bags (48 x 58cm) are heavy duty, high impact biohazard bags made of a 0.055mm thick high molecular weight polypropylene blend.

They are autoclavable to 138°C, saving autoclave time.

Each autoclave bag is printed with the standard biohazard warning symbol and precautionary procedures in four languages: English, Spanish, French and German.

They also feature a sterilization indicator patch that darkens when exposed to steam sterilization. Do not close the bags hermetically when placing them in the autoclave.

Specifications:

Dimensions (cm)	48 x 58
Volume (L)	24.4
Material	Polypropylene
Quantity	200



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com