

Deltalab Autoclave Resistant Red Bags (65 x 91 cm)

Product Code: 200321

Brand: Deltalab

Deltalab's Autoclave Resistant Red Bags (65 x 91 cm) are heavy duty, high impact biohazard bags made of a 0.055mm thick high molecular weight polypropylene blend.

They are autoclavable to 138°C, saving autoclave time.

Each autoclave bag is printed with the standard biohazard warning symbol and precautionary procedures in four languages: English, Spanish, French and German.

They also feature a sterilization indicator patch that darkens when exposed to steam sterilization. Do not close the bags hermetically when placing them in the autoclave.

Specifications:

Dimensions (cm)	65 x 91
Volume (L)	75.1
Material	Polypropylene
Quantity	200



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com