

## Deltalab Autoclave Resistant Red Bags (78 x 96cm)

**Product Code:** 200322

**Brand:** Deltalab

Deltalab's Autoclave Resistant Red Bags (78 x 96cm) are heavy duty, high impact biohazard bags made of a 0.055mm thick high molecular weight polypropylene blend.

They are autoclavable to 138°C, saving autoclave time.

Each autoclave bag is printed with the standard biohazard warning symbol and precautionary procedures in four languages: English, Spanish, French and German.

They also feature a sterilization indicator patch that darkens when exposed to steam sterilization. Do not close the bags hermetically when placing them in the autoclave.

### Specifications:

<b>Dimensions (cm)</b>	78 x 96
<b>Volume (L)</b>	121
<b>Material</b>	Polypropylene
<b>Quantity</b>	200



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)